



HOW TO CONFIGURE YOUR SYSTEM

CoolLED pE-2 excitation system

4 wavelength LED system with exchangeable dichroic mirror

1 Select up to 4 LED wavelengths - λ^*




Collimator fits directly to microscope's epi-fluorescent port.

2 CoolLED will provide pE-2 Collimator with a dichroic mirror optimised for LED wavelengths selected. Additional dichroic mirrors can be exchanged when alternative LAMs are installed.

3 Specify microscope. CoolLED provides adaptors to all current models.

4 For additional LAMs and accessories www.CoolLED.com/Products/



Product code	Description	Quantity
244-3750	pE-2 Main Unit (Drives 4 λ)	1
244-1500	Control Pod	1
244-4300	pE-2 Exchangeable Collimator	1
244-4312	pE-2 Exchangeable dichroic mirror	1
244-1904	Light Guide (1m)	2
244-2200	Microscope Adaptor (Specify)	1
244-5100- λ_1 - λ_2	LAM @ λ_1 + λ_2 nm	1
244-5100- λ_3 - λ_4	LAM @ λ_3 + λ_4 nm	1